

Sterilisation

- Autoclav (135°C)
- Heißluft (180°C)
- Thermo-Desinfektor



Siehe RÖDER Pflegeanweisung

Sterilization



- Autoclave (135°C / 266°F)
- Hot Air (180°C / 355°F)
- Thermo-Disinfection

Please see RÖDER cleaning instructions



Ausführungen

- Glas: plan 
- hohl 
- Silikoneinlage zur Abdichtung
- Stift: ss = simple stem
cs = cone socket
- rostfreier Stahl



Performances

- Glass: plane 
- magnifying 
- Silicone insert for sealing up
- Stem: ss = simple stem
cs = cone socket
- stainless steel

Verspiegelung

- Rückseitige
Verspiegelung 
- Oberflächen-
Verspiegelung 

Reflection

- Backside reflection 
- Front surface
reflection 

Verpackungseinheit

- Mundspiegel:
1 Packung enthält 12 Mundspiegel
- Spiegel für HNO:
Einzelverpackung

Packaging

- Mouth Mirrors:
1 package contains 12 mouth mirrors
- Laryngeal mirrors:
Single packaging

